

A Hoof Care Instruments

Page 1 - 8

B Farrier Tools

Page 9 - 12

C Castration Instruments

Page 13 -

Wide Hoof Knives Double Blade

Perfectly Sharpened Stainless Steel Blades With Matt Finish And Precise Quality Wooden Handle With Nice Finish Provides Better Grip.



Ref: 92.12.12



Ref: 92.12.13



Ref: 92.12.14



Ref: 92.12.15

Ref.	Type	Details
92.12.12	Right	Matt Finish, Standard Wood Handle
92.12.13	Left	Matt Finish, Standard Wood Handle
92.12.14	Right	Matt Finish, Special Wood Handle
92.12.15	Left	Matt Finish, Special Wood Handle

Single Edge Hoof Knives

Stainless Steel Blades With Matt Finish Are Sharpened With Extra Care So Provide Super Quality Cutting While Nicely Designed Handle From Special Quality Wood With Extra Finish Provides Excellent Grip.



Ref: 92.12.16



Ref: 92.12.17

Ref.	Type	Details
92.12.16	Right	Matt Finish, Fancy Wood Handle
92.12.17	Left	Matt Finish, Fancy Wood Handle

All Hoof Knives Models Are Available and can be customized upon your specifications i.e. Blades can be coated OR Finishing and handle design can be changed etc.

Customization Charges May Apply

Narrow Hoof Knives

Extra Narrow Highly Sharpened Blades Have Very Good Result In Cutting And Nicely Finished Fine Quality Wooden Handle Is Bet For Better Grip.



Ref: 92.12.18



Ref: 92.12.19

Ref.	Type	Details
92.12.18	Right	Matt Finish, Special Wood Handle
92.12.19	Left	Matt Finish, Special Wood Handle

Broad Hoof Knives

Wide Stainless Steel blades are perfectly sharpened to cut comfortably while better grip is assured due to very professionally prepared wooden handle.



Ref: 92.12.20



Ref: 92.12.21

Ref.	Type	Details
92.12.20	Right	Matt Finish, Fancy Wood Handle
92.12.21	Left	Matt Finish, Fancy Wood Handle

Narrow Hoof Knives Black Handle

Razor Edge Sharpen Steel Polished Blades Have Extra Hardness To Provide Maximum Sharpness. In Addition Prime Quality Black Rosewood Handle With Extra Thick Pins Gives More Durability And Powerful Grip.



Ref: 92.12.22



Ref: 92.12.23

Ref.	Type	Details
92.12.22	Right	Polished Razor Blade , Rosewood Handle
92.12.23	Left	Polished Razor Blade , Rosewood Handle

Narrow Hoof Knives

Extra Narrow Highly Sharpened Blades Have Very Good Result In Cutting And Nicely Finished Fine Quality Anti Slip **Plastic Handle** Is Best For Better Grip.

Ref: 92.12.24



Ref: 92.12.25



Ref.	Type	Details
92.12.24	Right	Matt Finish, Plastic Handle
92.12.25	Left	Matt Finish, Plastic Handle

Broad Hoof Knives

Wide Stainless Steel blades are perfectly sharpened to cut comfortably while better grip is assured due to very professionally prepared **Plastic handle**.

Ref: 92.12.28



Ref: 92.12.29



Ref.	Type	Details
92.12.28	Right	Matt Finish, Plastic Handle
92.12.29	Left	Matt Finish, Plastic Handle

Wide Hoof Knives Double Blade

Perfectly Sharpened Stainless Steel Blades With Matt Finish And Precise Quality **Plastic Handle** With Nice Finish Provides Better Grip.

Ref: 92.12.26



Ref: 92.12.27



Ref.	Type	Details
92.12.26	Right	Matt Finish, Plastic Handle
92.12.27	Left	Matt Finish, Plastic Handle

Cottam Hoof Knives Double Blade

Perfectly Sharpened Stainless Steel Blades With Polish Finish And Precise Quality Wooden Handle With Nice Finish Provides Better Grip.



Ref: 92.12.31



Ref: 92.12.32



Ref.	Type	Details
92.12.31	Right	Polish Finish, Multicolor Wooden Handle
92.12.32	Left	Polish Finish, Multicolor Wooden Handle

Finger Set Plastic Handel
Knife Image will be Here
92.12.30

Ascot Hoof Knives **Multicolor Handle**
Razor Edge Sharpen Blades Polished Have Extra Hardness To Provide Maximum Sharpness. In Addition Prime Quality **Multicolor Wooden Handle** With Extra Thick Pins Gives More Durability And Powerful Grip.

	Ref.	Details
1	92.12.33	Razor Sharpness , Polish Finish, Extra Fine Quality Wooden Handle
2	92.12.34	Standard Sharpness , Matt Finish, Standard Quality Wooden Handle

Frost Hoof Knife Black Coated

High Temper Stainless Steel Blades With Black Coating Are Very Nicely Sharpened So Provide Better Cutting While White Wood Handle With Extra Smoothness Provide Very Comfortable And Strong Grip.



Ref: 92.12.35

Ref.	Details
92.12.35	Razor Sharpness, Black Coated, Extra Fine Quality Wooden Handle

Frost Hoof Knife

High Temper Stainless Steel Blade is Very Nicely Sharpened So Provide Better Cutting While White Wood Handle With Extra Smoothness Provide Very Comfortable And Strong Grip.



Ref.	Details
1 92.12.36	Razor Sharpness, Polish Finish, Extra Fine Quality Wooden Handle
2 92.12.37	Standard Sharpness, Matt Finish, Standard Quality Wooden Handle

Small Hoof Knife

Steel Blades With Matt Finish Are Very Nicely Sharpened And Provide Good Quality And Easy Cutting While Wood Handle With Very Fine Finishing And Shiny Coating Provide Smoothness And Strong Grip.



Ref: 92.12.38

Ref.	Details
92.12.38	Matt Finish, Fine Quality Wooden Handle

Straight Hoof Knife Double Blade

Steel Blades With Matt Finish Are Very Nicely Sharpened And Provide Good Quality And Easy Cutting While Wood Handle With Very Fine Finishing And Shiny Coating Provide Smoothness And Strong Grip.



Ref: 92.12.39

Ref.	Details
92.12.39	Matt Finish, Fine Quality Wooden Handle

Bot Hoof Knife With Serrated Blade

Extra Hardness Of Stainless Steel Blades With Serrations Provide Supper Cutting And Wooden Handle Is Precisely Designed For Better Grip.



Ref: 92.12.40

Ref.	Details
92.12.40	Serrated Blade, Matt Finish, Fine Quality Wooden Handle

Cobra Hoof Knife

Cobra Shaped Matt Finished Blade Is Made From Good Quality Stainless So Have Wonderful Result In Cutting Also Wooden Handle With Special Pattern Is Providing Non Slip Grip By Using Less Force.



Ref: 92.12.41

Ref.	Details
92.12.41	Best Quality Sharpness with Matt Finish, Fine Quality Wooden Handle

Beak Hoof Knife

Supper Cutting Blade Having Two Beaks And Matt Finish Along With Wooden Handle Providing Strong Grip Makes Complex Cutting Of Hoof's Edges Very Easy Quick



Ref: 92.12.42

Ref.	Details
92.12.42	Two Beaks and Blade with Extra Sharpness and Matt Finish, Fine Quality Wooden Handle

Blunt Tip Hoof Knife

Stainless Steel blade with Wooden Handle and cutting adjuster removable.



Ref: 92.12.43

Ref.	Details
92.12.43	Stainless Steel Blade With Cutting Adjuster and Fine Quality Wooden Handle

Wide Loop Hoof Knife

This Wide Loop Hoof Knife Is Made From Superiors Class Stainless Steel And Is Both Side Sharpen To Get More Sharpness So Provide Better Result And Wooden Handle With A Special Coating Provide Very Nice Grip.

Ref: 92.12.44



Ref: 92.12.45



Ref.	Details
92.12.44	Matt Finish Stainless Steel Blade, Polished Wooden Handle
92.12.45	Mirror Finish Stainless Steel Blade, Multicolor Polished Wooden Handle

Medium Loop Hoof Knife

This Loop Hoof Knife Is Made From Superiors Class Stainless Steel And Is Specially Sharpened To Get More Sharpness So Provide Better Result And Wooden Handle With A Special Coating Provide Very Nice Grip.

Ref: 92.12.46



Ref.	Details
92.12.46	Matt Finish Stainless Steel Blade, Polished Wooden Handle

Small Loop Hoof Knife

This Loop Hoof Knife Is Made From Superiors Class Stainless Steel And Is Specially Sharpened To Get More Sharpness So Provide Better Result And Wooden Handle With A Special Coating Provide Very Nice Grip.

Ref: 92.12.47



Ref.	Details
92.12.47	Matt Finish Stainless Steel Blade, Polished Wooden Handle

Wide Loop Hoof Knife With Beak

This Wide Loop Hoof Knife Has A Beak On His Blade That Is Professionally Razor Sharpened, In Addition Teflon Handle Keep It Light Weight While The Brass Pins Provide More Durability And Its Extra Smoothness And Handle Shape Is Very Helpful To Provide Better Nonslip Grip.

Ref: 92.12.48



Ref.	Details
92.12.48	Razor Sharpness, Mirror Finish, Heavy Duty Pattern Teflon Handle With Brass Pins

Swiss Hoof Knife Standard Quality

Well Sharpened Loop Of This Swiss Hoof Knife Provide Nice Result In Cutting And Cleaning Of Hooves, Nicely Designed Wooden Handle Make It Easy To Use With

Ref: 92.12.49



Ref.	Details
92.12.49	Matt Finish Stainless Steel Blade, Polished Wooden Handle

Swiss Hoof Knife Premium Quality

Well Sharpened Loop Of This Swiss Hoof Knife Provide Nice Result In Cutting And Cleaning Of Hooves, Nicely Designed Wooden Handle Make It Easy To Use With Healthy Grip.

1 2



Ref.	Details
1 92.12.50	Matt Finish Stainless Steel Blade, Polished Wooden Handle
2 92.12.51	Matt Finish Stainless Steel Blade, Plastic Handle

All Loop Knives Models Are Available and can be customized upon your specifications i.e. Blades can be coated OR Finishing and handle design can be changed etc.

Customization Charges May Apply

Folding Hoof Knife With Hoof Pick
Stainless Steel Hoof Knife Blade And Hoof Pick Fold Able In One Very Durable Stainless Steel Handle. Easy To Open And Fold Back While Lock Able During Use.



Ref: 92.12.52

Ref.	Details
92.12.52	Lock Able During Use, Matt Finish Hoof Knife Blade And Hoof Pick With Matt Finish Handle.

Folding Knife For Hoof, 2 In 1
Stainless Steel Cutting Knife Blade And Pick Fold Able In One Very Durable Stainless Steel Handle. Easy To Open And Fold Back While Lock Able During Use.



Ref: 92.12.53

Ref.	Details
92.12.53	Lock Able During Use, Matt Finish Knife Blade And Pick With Matt Finish Handle.

Horseman's Multipurpose Knife
Five Blades Are Made Of Tempered Stainless Steel And Snap Securely In To The Rugged Wood And Steel Handle.



Ref: 92.12.54

Ref.	Details
92.12.54	Stainless Steel Blades and Handle With Wooden Shell. Matt Finish Blades. Knife Includes:- Leather awl, knife blade, hoof knife, hoof pick and thinning blade.

Hoof Knives Set of Seven with Universal Interchangeable Handle
Six blades are made of tempered stainless steel and snap securely in to the rugged wood and steel handle.

Set Includes:- Handle, Hoof Searcher, Sage Knife, Oval Hoof Knife Right Hand Wide, Right Hand Narrow, And Left Hand Wide Knife. **Leather carrying case.**



Ref.	Details
92.12.55	Usually Sharp Stainless Steel Blades, Mirror Finish And Stainless Steel Handel With Wooden Shell.
92.12.56	Razor Sharp Stainless Steel Blades, Mirror Finish And Stainless Steel Handel With Wooden Shell.

1
2

Hoof Buffer

It Has Several Different Uses And The Wide Blade Is Beveled To A Sharp Edge. The Over-all Length Is About 6" And The Wide Blade Is About 1 1/2" Wide.



Ref: 92.12.57

Ref: 92.12.58

Ref.	Details
92.12.57	100% Stainless Steel, Matt Finish, Overall Length is About 6 Inches, Wide Blade is About 1 1/2" Wide
92.12.58	100% Stainless Steel, Matt Finish, Overall Length is About 7 Inches, Wide Blade is About 1 3/4" Wide

Hoof Care Instruments

Hoof Picks, Scraper, Pincer, Cutter and Scissors

White Toe Nail Scissors

Made Of Fine Quality Stainless Steel, Dull/matt Finish.

Ref: 92.12.59



Ref.	Details
92.12.59	Stainless Steel, Matt Finish Size: 4.75 Inches

Cat Claw Scissors

Made Of Fine Quality Stainless Steel, Dull/matt Finish.

Ref: 92.12.60




Ref.	Details
92.12.60	Stainless Steel, Matt Finish Size: 4.75 Inches

Hoof Pick Small

Stainless Steel With Rubber Grip On Handle.

Ref: 92.12.61



Ref.	Details
92.12.61	Stainless Steel with Rubber Coated Grip on Handle All Colors are Available

Hoof Scraper

Razor Sharpened, Stainless Steel with Rubber Dipped Grip on Handle.

Ref: 92.12.63



Ref.	Details
92.12.63	Razor Sharpened, Stainless Steel, Rubber Grip on Handle

Heavy Duty Cutter

Carbon Steel Made With Black And Red Coating And Rubber Grip On Handle, Professionally Sharpened And Provide Maximum Result In Cutting, Easy To Use

Ref: 92.12.64



Ref.	Details
92.12.64	Carbon Steel, Black Coated Blades, Red Coated Body and Rubber Grip on Handle

Hoof Pick Large

Fine Quality Stainless Steel, Dull Finish. Plastic Molded Handle.

Ref: 92.12.62



Ref.	Details
92.12.62	Stainless Steel with Plastic Molded Handle

Hoof Pincer

High Temper Stainless Steel Made, Very Sharpened Blades That Make Cutting Easy And Long Time.

Ref: 92.12.65



Ref.	Details
92.12.65	Stainless Steel, Professionally Sharpened, Easy Cutting

Foot Root Hoof Shear

Die Caste Forged Stainless Steel Made With Matt Finish Blades Are Very Resistant And Sharp. Spring Is With Extra Durability. Shear Is Lock Able With Leather Strip Lock.

Ref: 92.12.66



Ref: 92.12.67



Ref.	Size	Details
92.12.66	8 Inch	Stainless Steel with Matt Finish & Lock
92.12.67	8 Inch	Serrated Matt Finish Blade & Black Coated Handle With Lock

Foot Root Hoof Shear Heavy Duty

High Tempered Stainless Steel Drop Forged Blade and Die Casted Handle are combined with rivets. Mirror finish provides very soft but strong grip. Handle can be locked.

Ref: 92.12.68



Ref: 92.12.69



Ref.	Size	Details
92.12.68	8 Inch	Stainless Steel with Mirror Finish & Lock
92.12.69	8 Inch	Serrated Matt Finish Blade & Green Coated Handle with Lock



92.12-70 Blister Packed
92.12-71 Seal Packed

Heavy Duty Hoof Trimmers

High Quality Stainless Steel with Maximum Hardness & Blades with Extra Sharpness to Cut Easily. Rubber Grips on handles for Strong and Comfortable Holding



Ref: 92.12.72



Ref: 92.12.73

Ref.	Size	Details
92.12.72	18 Inch	Standard Quality, Carbon Steel with Coating and Rubber Grip on Handle
92.12.73	18 Inch	Prime Quality, Carbon Steel with Green Coating and Rubber Grip on Handle

Hoof Testers

Fine Quality High Temper Stainless Steel with Dull Finishing Look very impressive as well is very best to detect hoof injuries easily.



Ref: 92.12.74



Ref: 92.12.75



Ref: 92.12.76



Ref: 92.12.77



Ref: 92.12.78



Ref: 92.12.79

Ref.	Size	Details
92.12.74	13" (inch)	Matt Finish, Stainless Steel
92.12.75	15" (inch)	Adjustable, Matt Finish, Stainless Steel
92.12.76	15" (inch)	Matt Finish, Stainless Steel
92.12.77	16" (inch)	Matt Finish, Stainless Steel
92.12.78	16" (inch)	Matt Finish, Stainless Steel
92.12.79	19" (inch)	Heavy Duty, Matt Finish, Stainless Steel

Hoof Nippers

Made from Chrome Vanadium Steel by Drop Forging Method and are precision machined to obtain the correct bevel to the cutting edge. Extra Hard Work on the alignment and sharpening of blades results accurate trimming and long life to the blades while excellent design are best for professional farriers to use any time with maximum precision and minimum fatigue.

Coatings And Packing
PVD (titanium), Jet Black And Other Types Of Coatings And Special Bespoke Packing Is Available Upon Your Requirements At Additional Charges

Ref.	Size	Details
1	92.13.01	12" Matt Finish, Chrome Vanadium Steel
2	92.13.02	14" Matt Finish, Chrome Vanadium Steel
3	92.13.03	15" Matt Finish, Chrome Vanadium Steel
4	92.13.04	12" Matt Finish, Chrome Vanadium Steel
5	92.13.05	14" Matt Finish, Chrome Vanadium Steel
6	92.13.06	15" Matt Finish, Chrome Vanadium Steel
7	92.13.07	12" Matt Finish, Chrome Vanadium Steel
8	92.13.08	14" Matt Finish, Chrome Vanadium Steel
9	92.13.09	15" Matt Finish, Chrome Vanadium Steel
10	92.13.10	12" Jet Black, Chrome Vanadium Steel
11	92.13.11	14" Jet Black, Chrome Vanadium Steel
12	92.13.12	15" Jet Black, Chrome Vanadium Steel
13	92.13.13	12" Jet Black, Chrome Vanadium Steel
14	92.13.14	14" Jet Black, Chrome Vanadium Steel
15	92.13.15	15" Jet Black, Chrome Vanadium Steel



Pull Off

Made from Chrome Vanadium Steel by Drop Forging Method and are precision machined to obtain the correct bevel to the edge to make shoe removal easier and also cut horseshoe nails without damage to the blades.

Ref.	Size	Details
1	92.13.16	15" Side Teeth, Matt Finish, Chrome Vanadium Steel
2	92.13.17	12" Heavy Head Matt Finish, Chrome Vanadium Steel
3	92.13.18	14" Heavy Head Matt Finish, Chrome Vanadium Steel
4	92.13.19	12" Normal Head Matt Finish, Chrome Vanadium Steel
5	92.13.20	14" Normal Head Matt Finish, Chrome Vanadium Steel
6	92.13.21	15" Side Teeth, Matt Finish, Chrome Vanadium Steel



Customization
All Products Can Be Customized As Per Your Requirement Upon You Specifications, Drawings, Samples, Images Or Reference
Customization Charges May Apply

Clinchers

Made From High Quality Tools Steel By Drop Forging Method And Are Precision Machined With Hard Work On The Alignment To Increase The Functionality And Performance.

Teeth Are Heat Treated To Hold Their Bite Longer Under Heavy Use.

All Clinchers Are Easy To Use On All Horses.

Ref: 92.13.22



Ref: 92.13.23



Ref: 92.13.24



Ref: 92.13.25



Ref: 92.13.26



Ref: 92.13.27



Ref.	Size	Details
92.13.22	12 ½"	Clincher, Matt Finish, High Quality Tools Steel
92.13.23	10"	High Nail Clincher, Matt Finish, Quality Tools Steel
92.13.24	10"	High Nail Clincher, Matt Finish, Quality Tools Steel
92.13.25	13"	Master Clincher, Matt Finish, Quality Tools Steel
92.13.26	12 ½"	Clincher, Jet Black Finish, High Quality Tools Steel
92.13.27	12"	Clincher, Jet Black With Rubber Grip on Handle


Nail Puller

As Our Other Farrier Tools Our Nail Puller Are Also Drop Forged From Very High Quality Tools Steel. The Jaws Are Hardened To Give The Highest Durability Are Designed To Lock Over The Nail Head And Along Its Taper To Easily Remove A Single Nail From A Creased Swedged, Or Fullered Shoe Very Easily.

Ref: 92.13.28



Ref: 92.13.29



Ref.	Size	Details
92.13.28	12"	Nail Puller, Matt Finish, High Quality Tools Steel
92.13.29	12"	Nail Puller, Matt Finish, High Quality Tools Steel

Shoeing Tong

You Can Hold Draft Shoes With Ease. The Hinge Works Beautifully, And You Can Use Them Either With The Prong End, Or With The Other End And You Have A Good Hold On A Hot Shoe.

This Shoeing Tong Is Made From 1045 Steel Because It Can Withstand Repeated Heating And Cooling. It Is Hand Forged At The Anvil And Is Fit To The Horseshoe, Allowing A One Handed Grip. This Tool Is Designed To Be Used On Either End.

Ref: 92.13.30



Ref.	Size	Details
92.13.30	15"	Shoeing Tong, Matt Finish, High Quality Tools Steel

Coatings And Packing
PVD (titanium), Jet Black And Other Types Of Coatings And Special Bespoke Packing Is Available Upon Your Requirements At Additional Charges

Heavy Duty Hammers CNC Crafted

Manufactured From High Quality Alloy Steel. 12" Long Ergonomic Hardwood Well Balanced Handle.

This Tool Boasts One Rounded Face And One Flat Face With Rounded Edges To Help Eliminate Hammer Marks When The Steel Is Not Hit Exactly Flat.

Further These Hammers Are Tempered Back To Be Slightly Softer Through The Bodies, And Then Hardened Only On The Striking Ends. This 3 Step Process Takes More Time, But Results In Less Vibration, And Shock, With Each Blow Of The Hammer.



Ref: 92.13.31



Ref: 92.13.32



Ref: 92.13.33



Ref: 92.13.34



Ref: 92.13.35

Ref.	Weight	Details
92.13.31	1000 Grams	Matt Finish, 12" Long Ergonomic Hardwood Handle
92.13.32	1000 Grams	Matt Finish, 12" Long Ergonomic Hardwood Handle
92.13.33	1000 Grams	Matt Finish, 12" Long Ergonomic Hardwood Handle
92.13.34	1000 Grams	Matt Finish, 12" Long Ergonomic Hardwood Handle
92.13.35	1000 Grams	Matt Finish, 12" Long Ergonomic Hardwood Handle

Hoof Hammers

Made from high grade forged steel. 12" long ergonomic hickory handle and well balanced. The jaws are designed to ring off nails with ease. Hammer can also be used to remove nails that have been well driven.



Ref: 92.13.36



Ref: 92.13.37



Ref: 92.13.38



Ref: 92.13.39

Ref.	Weight	Details
92.13.36	350 Grams	Matt Finish, 12" Long Ergonomic Hickory Handle
92.13.37	350 Grams	Matt Finish, 12" Long Ergonomic Hickory Handle

Ref.	Weight	Details
92.13.38	350 Grams	Matt Finish, 12" Long Ergonomic Hickory Handle
92.13.39	350 Grams	Matt Finish, 12" Long Ergonomic Hickory Handle

Hoof Hammers With Springs

Made From High Grade Forged Steel. 12" Long Ergonomic Hickory Handle And Well Balanced. The Jaws Are Designed To Ring Off Nails With Ease. Supports Keeps Handle Tight And Provide Little Extra "whip". Hammer Can Also Be Used To Remove Nails That Have Been Well Driven.



Ref: 92.13.40



Ref: 92.13.41

Ref.	Weight	Details
92.13.40	350 Grams	Matt Finish, 12" Long Ergonomic Hickory Handle
92.13.41	350 Grams	Matt Finish, 12" Long Ergonomic Hickory Handle

Hoof Rasps Without Tang

Professionally Made with 4 sides with fine rasp cut

Ref: 92.13.42



Ref: 92.13.43



Ref.	Details
92.13.42	Prime Quality Rasp with Four Side Fine Quality Cuts
92.13.43	Prime Quality Rasp with Four Side Fine Quality Cuts

Curved Hoof Rasps Without Tang

Two-sided With Fine And Rough Files

Ref: 92.13.44




Ref.	Details
92.13.44	One Side Fine And One Side Rough Cuts

High End Hoof Rasps

Durability Up To 3x. One Side Extra Sharp And Deep Rasp Teeth For With Excellent Performance While Other Revolutionary And Distinctive File Side For A Fine Surface. Comfortable Smooth Running

Ref: 92.13.45



Ref.	Details
92.13.45	2 Sided, One Sharp and Deep, Other File Type

Hoof Rasps Heavy Duty

Heavy Duty Hoof Rasp Is Designed Specifically For Smoothing And Filing Hooves. Coarse Teeth On One Side Fine Teeth On Other. 14" Long.

Ref: 92.13.46



Ref.	Details
92.13.46	2 Sided, One Deep Sharp, Other Regular Sharp

Hoof Rasps Large

Professional Quality Rasp Has Rasp Teeth On One Side And File Teeth On The Other Side, Single Cut On Edges

Ref: 92.13.47



Ref.	Details
92.13.47	2 Sided, One Rasp Teeth, Other File Teeth

Hoof Rasps

Designed To Perform Well In Both Wet And Dry Conditions. Angled Tooth Face Resists Clogging. Rasp Teeth On One Side, File Teeth On The Other Side. Single Cut On Edges

Ref: 92.13.48



Ref.	Details
92.13.48	2 Sided, One Side Rasp and Other Side File Teeth

Wooden Handle For Hoof Rasps

Made From Special Quality Hard Wood Handle Fits Easily Onto Rasp

Ref: 92.13.49 **Large**



Small Ref: 92.13.50



Ref: 92.13.51 **Universal**



Leather Cover For Hoof Rasps

Fine Quality Handy Rasp Cover To Keep Rasps Clean And Sharp

Ref: 92.13.52

